



# KDZ Bruxer Aesthetic®

Highly Translucent Anterior Zirconia



*Keating Dental Lab has provided dental practices with the highest quality restorations, paired with reliable, friendly service since 2002.*

*As a family-run full-service dental lab, we treat our dentists like we would our own. Our goal is to provide doctors with the expertise of a national lab, and the convenience of a local laboratory partner.*



## THE MOST LIFELIKE ZIRCONIA CROWN AND BRIDGE OPTION

Bring the reliable strength of full-contour zirconia into the cosmetic zone with KDZ Bruxer Aesthetic® crowns and bridges. Keating Dental Lab's revolutionary new zirconia formulation provides the ultimate in translucency, rivaling IPS e.max® crowns, while its 650 MPa flexural strength is nearly twice that of lithium disilicate.

### INDICATIONS FOR KDZ BRUXER AESTHETIC®:

- Ideal for anterior and premolar cases
- Single unit crowns and bridges with one pontic
- 0.8mm of reduction required

*CAD/CAM fabricated to match all 16 VITA® Classic shades with the industry-leading precision of our digital workflow, KDZ Bruxer Aesthetic® is the next-generation zirconia perfect for demanding patients.*

