

KDZ Bruxer®

Full Contour Zirconia Crowns



Keating Dental Lab has provided dental practices with the highest quality restorations, paired with reliable, friendly service since 2002.

As a family-owned and operated lab, we treat our dentists like we would our own. Our goal is to provide doctors with the expertise of a national lab, and the convenience of a family-run lab.



CERAMIC DENTAL IMPLANT CROWNS – KDZ BRUXER®

With zirconia formulated by Keating Dental Lab's research and development team to meet our demanding standards for durability and aesthetics, the KDZ Bruxer® crowns are the most esthetic monolithic zirconia restorations on the market. They are gentle on opposing enamel for bruxing patients, yet strong enough for long-span bridges with an average flexural strength of 1,100 MPa. The KDZ Bruxer® crowns are available in all 16 Classic Vita shades.

INDICATIONS FOR KDZ BRUXER:

- Posterior Crowns or Bridges
- Screw-Retained Crowns or Bridges
- Crowns or Bridges over Custom Abutments

An extreme strength of CAD/CAM monolithic Zirconia that is an esthetic alternative to the traditional Full Cast restoration. Three to six times the fracture toughness of typical zirconia. Can be conventionally cemented using non-expanding resin modified glass ionomer.

